Expt. No. 2 Page No. _ 2 Exporiment - 2 Estimation of Hardness of Water → Objective:- To determine total, permanent & temporary handness of given water sample by complexometric titration using EDTA sol > Psinciple: - The presence of ca+2 & un water introduces hardness. Hard water consumes lot of soop before forming from the sound water is used in a tailer, it hereto in crossive scaling & sludging which result in loss of efficiency and I cate & Mgt2 wons are present as a carbonate salts, they can be easily removed by boiling, decomposes bicarcomete to in aluble cartonale salts this is called temporary When Cat2 & Mg+2 are hardness. present in any other form it is called berman ent hardness To determine hardness we pedorm titration. EBT is used as indicator and ditration is performed with EDTA Refer to reaction (i Teacher's Signature ___

Date 10/2/2025

W	nen	EBT (ů oc	td ad	to hand	
- AUG	NUMBER	011		10440	least a	1 1 1
LAGA		The last of	(I) no	FT	A Comment	66 6 6
-414		-	COX (094 Y D)	A Attack 1	14 40	- 2.4
- 11	cu	MOIN	4 0	0	WOLLENS .	12 EAST
ubo	come	s fore	se .	- Hence	end boil	it of
th	1	titicot	ion c	is ide	end pointified to	color.
Chi	ange	Re	fer .	to 4	teaction (ii)	
of	0	Adution	nea	ely	l to	pH - 10
Ob.	seewa	tions				
Ob.	y H	andness	Burette	Reading	Val of EDTA	Ava Burette
Ob.	y H	tions	Burette	Reading		Ava Burette
Ob.	Vol Har	asidness a waron(mi)	Buxette Initial 37.5	Reading at Final) val of EDTA Consumed (m1)	Avg Burette Reading
Ob.	You Har	ardness of a waror(m)	Buxette Initial 37.5	Readings i	Val of EDTA Consumed (m1) 12-1 12-7	Avg Burette Reading
Ob.	You Har	ardness of a waror(m)	Buxette Initial 37.5	Reading at Final	Val of EDTA Consumed (m1) 12-1 12-7	Avg Burette Reading